


MEIER-AUGENSTEIN, W. and SCHIMMELMANN, A. 2019. A guide for proper utilisation of stable isotope reference materials. *Isotopes in environmental and health studies* [online], 55(2), pages 113-128. Available from: <https://doi.org/10.1080/10256016.2018.1538137>

A guide for proper utilisation of stable isotope reference materials.

MEIER-AUGENSTEIN, W. and SCHIMMELMANN, A.

2019

This is an Accepted Manuscript of an article published by Taylor & Francis in *Isotopes in environmental and health studies* on 08-11-2018, available online:
<http://www.tandfonline.com/10.1080/10256016.2018.1538137>.



A guide for proper utilization of Stable Isotope Reference Materials

Journal:	<i>Isotopes in Environmental & Health Studies</i>
Manuscript ID	Draft
Manuscript Type:	Original Article
Date Submitted by the Author:	n/a
Complete List of Authors:	Meier-Augenstein, Wolfram; Robert Gordon University, School of Pharmacy & Life Science Schimmelmann, Arndt; University of Indiana, Department of Geological Sciences
Keywords:	isotope abundance, Isotope measurements, methods and equipment, reference material, scale anchor, stretch factor, scale compression, scale normalisation, delta scale reference point

SCHOLARONE™
Manuscripts

A guide for proper utilization of Stable Isotope Reference Materials

Wolfram Meier-Augenstein¹ and Arndt Schimmelmann²

¹*School of Pharmacy and Life Sciences, Robert Gordon University, The Sir Ian Wood Building, Garthdee Road, Aberdeen, AB10 7GJ, UK; w.meier-augenstein@rgu.ac.uk*

²*Department of Earth and Atmospheric Sciences, Indiana University, 1001 East 10th Street, Bloomington, Indiana 47405-1405, United States; aschimme@indiana.edu*

A guide for proper utilization of Stable Isotope Reference Materials

Abstract

Many scientific publications about stable isotope ratios suffer from flawed practices regarding calibration and normalisation of raw δ values in conjunction with prescribed δ values of reference materials. Violations of the Identical Treatment principle with regards to samples and standards (i.e. reference materials) and lack of adherence to SI-mandated and IUPAC-recommended nomenclature exacerbate the widespread problem of lackadaisical analytical practice and reporting. Science is supposed to strive for exactness, whereas ambiguity and jargon confound interdisciplinary communication. This contribution aims to expose typical misconceptions and avoidable errors and offers guidance toward reproducible generation of isotope data, isotopic scale normalisation, and proper data reporting. We offer a comprehensive overview of sources of light stable isotope reference materials to best match sample matrices encountered by stable isotope practitioners with chemically similar reference materials.

Keywords: calibration; isotope abundance; isotope abundance calibration; reference material; scale anchor; scale compression; scale normalisation; stable isotope; stretch factor

1. Introduction

The exponential growth of global scientific output fosters compartmentalisation into sub-disciplines with specialized jargons and cryptic conventions. Practitioners of stable isotope analyses are spread over many poorly connected scientific fields from archaeology over ecology and forensics to geochemistry to name but a few [1-5]. Fruitful interdisciplinary communication mandates that scientific results should be measured reproducibly and reported unambiguously based on internationally accepted scientific units and nomenclature [6, 7]. This contribution aims to provide a guide towards generating and reporting stable isotope data with quality assurance and transparency. Adherence to such principles will facilitate publication of results as well as interdisciplinary understanding.

1
2
3 It is probably fair to say neither reviewers nor editors of scientific journals would
4 accept and publish manuscripts reporting quantitative data based on mass spectrometric (MS)
5 analysis if such data were not supported by a multi-point calibration. So, why should
6 manuscripts reporting isotope abundance data based on isotope ratio mass spectrometric
7 (IRMS) analysis be treated any differently? The answer of course is they should not. While
8 failure to calibrate an analytical method to compensate for method and/or instrument inherent
9 differences between measured and accepted (known to be true) values will present problems in
10 any subject of natural or life sciences, this is of particular concern in subject areas
11 impacting on public health and safety including consumer protection such as food
12 authenticity (or food forensics), authenticity of pharmaceutical drugs or forensic science.
13 Lack of traceability to internationally recognised scale anchors and lack of inter-laboratory
14 comparability of such non-calibrated data undermines not just confidence in the data but
15 confidence in the conclusions drawn. This would be particularly regrettable if such data were to
16 be presented and refuted in court [8]. Yet, many an example of articles reporting such data can be
17 found in journals dedicated to food authenticity or forensic science to this day.

18
19 This article therefore aims to dispel any misconceptions that may still exist
20 concerning isotope abundance calibration or isotopic scale normalisation of measured stable
21 isotope abundance data by explaining the concept of stable isotope reference materials and by
22 illustrating how to use them and how not to.

2. Why multi-point isotopic calibration is a must

23
24 Having referred to multi-point calibration being a prerequisite for quantitative MS analysis, let
25 us briefly explore the similarity or analogous principles between quantitative compound MS
26 analysis and compound isotope abundance analysis by IRMS. In quantitative MS analysis
27 calibration aims to determine the relation between compound amount or

1
2
3 concentration and detector response. Even though the detector response of MS instruments is
4
5 directly proportional to compound amount, performance variability of the different
6
7 instrument components results in uncertainty levels of > 10 % relative standard deviation
8
9 (RSD), too poor for reliable quantification. Repeatedly carried out calibration series can
10
11 therefore yield widely differing calibration curves even though individually each calibration
12
13 curve can be all but perfectly linear (Figure 1a). To overcome this problem quantitative MS
14
15 analysis relies on the use of internal standards (IS). Calibration curves are based on the
16
17 analysis of varying concentrations of the target compound(s) to which always the same
18
19 amount of internal standard(s) has been added. Calibration curves are constructed by plotting
20
21 abundance ratios of quantifying ion compound over quantifying ion IS (Figure 1b). This
22
23 approach is associated with uncertainty levels of < 1 % RSD since any factors ultimately
24
25 affecting detector response will affect target compound and IS measurement in the same way
26
27 and by building the ratio will cancel each other out.
28
29
30
31
32

33
34 The situation we encounter when analysing compounds for the abundance of the
35
36 heavier isotope of a given light element by IRMS is similar and yet somewhat different. One
37
38 could say calibration in quantitative MS and IRMS are two different sides of the same coin.
39
40 Isotope abundance analysis also aims to determine a quantitative answer, namely, the true
41
42 amount of a given isotope present in a compound or material. However, for a number of
43
44 reasons the analytical answer is not an absolute amount but a ratio of the heavier isotope
45
46 abundance over the lighter isotope abundance of the sample (S) relative to the known isotope
47
48 abundance ratio of a standard (STD).
49
50
51

$$\delta^h E_{S/STD} = \frac{R_S - R_{STD}}{R_{STD}} = \frac{R_S}{R_{STD}} - 1 \quad (1)$$

52
53
54
55
56 The result of this ratio of ratios calculation is the delta value (δ) of the heavier isotope
57
58 h of a given element E (Equation 1). Numerical results of this equation are typically quite
59
60

1
2
3 small, e.g. -0.02996 which is why for reasons of convenience they may be expressed as ‰
4
5 values where -0.02996 is written as -29.96×10^{-3} or -29.96 ‰. So, like quantitative MS
6
7 analysis of compound abundance IRMS analysis of stable isotope abundance yields a
8
9 quantitative answer yet with one subtle difference. Unlike quantitative MS analysis the result
10
11 of IRMS analysis can be either a negative or a positive number because measurement results
12
13 are expressed relative to an internationally agreed scale zero point. A negative number
14
15 therefore indicates the stable isotope abundance in a given compound is less than that of the
16
17 primary reference material defining the scale zero point. Conversely, a positive number
18
19 indicates the stable isotope abundance in a given compound is higher than that of the primary
20
21 reference material defining the scale zero point.
22
23
24

25
26 It is important to note, like ‰ values this is purely a mathematical convention but
27
28 should not be confused with an SI unit of measurement (SI = Système International a.k.a.
29
30 International System). Should one wish to avoid the ‰ notation (and thus the ‰ symbol)
31
32 being confused with an SI unit, particularly in graphs and tables, one can present relative
33
34 abundance data as $10^3 \times \delta$ values [9, 10]. In line with the established practice to attribute SI
35
36 supplementary units to values derived from ratios where units cancel each other (e.g. the
37
38 Reynolds number [Re]) the term "urey" (symbol Ur) has been proposed so that δ values
39
40 traditionally written as e.g. -29.96 ‰ can be written as -29.96 mUr [6].
41
42
43

44
45 While in theory any compound could be used as standard as long as its isotope
46
47 abundance with regards to ^1H were known or defined, it is obvious that δ values thus obtained
48
49 could not be compared between different laboratories unless all laboratories would use
50
51 subsamples of the same standard material. However, would one standard suffice to quantify
52
53 isotope abundance of a given element in a given compound or material? The answer is no of
54
55 course. We know ^{13}C abundance in chemically identical compounds like sugars can differ
56
57 widely depending by which photosynthetic pathway they were produced [11-13]. In other
58
59
60

1
2
3
4 words, in analogy to the need in quantitative MS analysis for concentration or abundance
5
6 calibration based on several but at least two standard preparations containing different
7
8 amounts of the target compound there is a similar need in quantitative isotope abundance
9
10 IRMS analysis for several but at least two standards of preferably identical chemical
11
12 composition but of different isotope abundance (i.e. isotopic specific quantity). For this
13
14 reason alone δ values obtained by comparative measurement to a single cylinder gas are
15
16 neither calibrated nor scale normalized. It could of course be argued multi-point calibration
17
18 should not be necessary considering we are measuring relative isotope abundance values.
19
20 However, such an argument fails to take into account a phenomenon particular to IRMS
21
22 analysis called scale compression [14, 15]. Scale compression is essentially the sum of all
23
24 mass discriminatory effects associated with sample gas transfer to the IRMS, sample gas
25
26 admission into the ion source of the IRMS and possibly processes inside the ion source itself.
27
28 The measurable consequence of scale compression is most noticeable and thus most
29
30 detrimental to accuracy in ^2H abundance analysis as illustrated by the examples presented in
31
32 Table 1. In that instance the difference between SLAP2's measured $\delta^2\text{H}$ value and its scale
33
34 normalized and thus isotopically calibrated $\delta^2\text{H}_{\text{VSMOW}}$ value was 31.8 ‰. Similarly, the
35
36 difference between the measured $\delta^2\text{H}$ value of NBS 22 and its scale calibrated $\delta^2\text{H}_{\text{VSMOW}}$
37
38 value of -117.0 ‰ (accepted $\delta^2\text{H}_{\text{VSMOW}}$ value of -117.2 ‰ [16]) was 8.5 ‰, representing a
39
40 relative difference of 7.25 % between measured and accepted ^2H abundance value. This
41
42 example alone illustrates the importance of using reference materials when carrying out
43
44 stable isotope abundance measurements. Equally, this example illustrates the importance of
45
46 reporting which international reference materials have been used as scale anchors for the
47
48 respective stable isotope scales and which accepted δ values have been employed [7]. There
49
50 are numerous instances where well-established international reference materials like NBS 22
51
52 or IAEA-CH-7 have been repeatedly isotopically re-evaluated over time.
53
54
55
56
57
58
59
60

3. Why internationally accepted reference materials are important

The inescapable conclusions from the foregoing are these. Internationally accepted reference materials (RMs) are required that satisfy two needs.

- (I) Ensuring relative isotope abundance measurements can be calibrated to internationally agreed methodology so calibrated δ values can be compared between laboratories on a like-for-like basis.
- (II) Ensuring relative isotope abundance measurements can be corrected for any scale distorting effects so calibrated delta values are anchored to internationally agreed isotope abundance scales as defined for a particular element.

In addition to delivering highly comparable data, reference materials meeting these two conditions will also put in place a good traceability system. There are other practical but equally important reasons why internationally accepted reference materials are needed in addition to the primary, scale defining reference materials [7] such as Vienna Standard Mean Ocean Water (VSMOW and VSMOW2) and Standard Light Antarctic Precipitation (SLAP and SLAP2) which scale anchor the ^2H as well as the ^{18}O δ scale [6, 10, 17, 18]. For one, for historical reasons the primary scale defining reference materials are inorganic compounds whose chemical elemental composition does not match that of organic compounds and materials. For another, limited stocks of primary reference materials place restrictions on how much and at which intervals they can be supplied to any one laboratory. To overcome these constraints, internationally accepted and distributed secondary reference materials have been made available "to bridge the materials and chemistry gap" [7]. Secondary reference materials are traceable to the scale defining primary reference materials but unlike the defined δ values of primary reference materials, specified δ values of secondary reference

materials are associated with an uncertainty envelope (Table 2). A list of distributors of reference materials is given in Table 3.

However, even secondary reference materials (Table 2) are not available in limitless supply, which means it is not practically feasible to use them directly as Equation 1 implies. They are therefore more efficiently used to scale normalise and thus calibrate "raw" δ values that have been measured against a working gas. It must be emphasized, measured δ values calculated by instrument software on the basis of a working gas peak are neither calibrated nor scale normalised even if the isotopic composition of the working gas is known. As illustrated by the example shown in Table 1, scale normalization requires at least two reference materials acting as scale anchors. Furthermore, in continuous-flow IRMS instruments working gas pulses are introduced in a gas stream separate to the sample gas stream and are not subject to the same physical and chemical processes as the samples. For these two reasons the use of working gas pulses does not meet with the Identical Treatment principle [7, 19, 20].

What does work however while meeting all requirements discussed above is using reference materials and working gases in combination (Figure 2; Table 4). This means both samples and reference materials serving as scale anchors are analysed under identical conditions and their measured δ values are expressed relative to a working gas. Measured δ values of reference materials 1 and 2 (RM1 and RM2) are subsequently compared to their accepted δ values on the relevant isotope reference scale. A linear regression line is fitted through the data points (RM1_{measured}, RM1_{accepted}) and (RM2_{measured}, RM2_{accepted}) to yield a scale normalization equation of the form

$$\delta^h E_{\text{scale}} = s \times \delta^h E_{\text{measured}} + b \quad (2)$$

$$\text{where } s = (\text{RM2}_{\text{accepted}} - \text{RM1}_{\text{accepted}}) / (\text{RM2}_{\text{measured}} - \text{RM1}_{\text{measured}})$$

$$\text{and } b = \text{RM1}_{\text{accepted}} - (s \times \text{RM1}_{\text{measured}}) = \text{RM2}_{\text{accepted}} - (s \times \text{RM2}_{\text{measured}})$$

1
2
3 Given what has been said thus far, it should be obvious that in principle the δ value of
4 the working gas does not have to be known and could be set arbitrarily to any value.
5
6 However, there are practical reasons why it is useful for the approximate δ value of the
7
8 working gas to be known. This is easily achieved by measuring pulses of the working gas
9
10 against a single reference material. The most important practical reason for doing so is for
11
12 raw δ values (i.e. based on working gas δ values) to serve as decision enabling acceptance
13
14 criteria. More often than not, sample material is in limited supply so a decision as to whether
15
16 a batch run analysis can proceed or should be aborted has to be made early on, e.g. based on
17
18 measured δ values of RM1 and thus before any samples (Table 4). In fact, standard operating
19
20 procedures for accredited analytical processes require defined acceptance criteria on which
21
22 basis a decision is made if an analytical batch sequence has to be stopped or is permitted to
23
24 proceed.
25
26
27
28
29
30
31
32
33
34

35 **4. Why and which secondary reference materials are appropriate scale anchors**

36
37 As discussed above, in terms of sample turn-around and efficient use of reference materials
38
39 working gases offer a convenient way to generate raw δ values indicative of a sample's actual
40
41 isotope abundance. Cylinder gases, however well their isotopic composition may be known
42
43 are not equivalent to the scale defining primary reference materials that anchor a particular
44
45 δ^{hE} value reference scale. Most scale defining primary reference materials are however
46
47 inorganic compounds or materials that are either not directly amenable or directly
48
49 comparable to stable isotope analysis of organic compounds and materials by Continuous
50
51 Flow - IRMS (CF-IRMS). For this reason, numerous secondary (often organic) reference
52
53 materials have been developed (Table 2). Like primary reference materials secondary
54
55 reference materials are internationally agreed-upon and are administered and distributed by
56
57 internationally
58
59
60

1
2
3 recognised organisations (Table 3). Their scale normalised δ values are based on statistically
4
5 valid results submitted by stable isotope laboratories that had participated in international
6
7 inter-laboratory exercises organised by the IAEA or other organisations. The purpose of
8
9 internationally distributed reference materials is to anchor δ^{hE} scales and to enable
10
11 comparable measurement results in stable isotope laboratories no matter their location in the
12
13 world. The majority of these secondary international reference materials are listed on the
14
15 website of the International Union of Pure and Applied Chemistry's (IUPAC) Commission on
16
17 Isotopic Abundances and Atomic Weights (CIAAW) at www.ciaaw.org. The importance of
18
19 international isotopic reference materials to produce robust, traceable and internationally
20
21 comparable results has been emphasized by IUPAC guidelines that have been widely
22
23 published [6, 7, 10].
24
25
26
27
28
29

- 30 • If a secondary international measurement standard defines the size of a δ^{hE} scale, such
31 as SLAP (SLAP2) water for $\delta^2\text{H}$ and $\delta^{18}\text{O}$ measurements, δ^{hE} values should be
32 normalised using both standards. The authors should state this clearly in their articles
33 and reports.
34
35
- 36 • Authors are also encouraged to analyse with their samples and report δ^{hE} values of
37 further internationally distributed reference materials as appropriate for the
38 measurement method concerned. Secondary internationally distributed isotopic
39 reference materials that are of a nature similar to those of the unknowns being
40 measured (sulphate, nitrate, cellulose, etc.) should be analysed. This has been called
41 the Identical Treatment principle [19] and minimises systematic errors by subjecting
42 sample unknowns and reference materials to exactly the same chemical and other
43 manipulation steps, including the transfer pathway to the ion source of the IRMS. In
44 this manner, measurement results can be adjusted in the future as analytical methods
45
46
47
48
49
50
51
52
53
54
55
56
57
58
59
60

1
2
3 improve and consensus values of internationally distributed isotopic reference
4 materials are amended.
5
6
7

8
9 Therefore it follows from the above preferably any two secondary reference materials
10 whose δ^hE values extend closely to the size of the relevant δ^hE scale can serve as appropriate
11 scale anchors. Even more preferably, such two secondary reference materials should be a
12 good matrix match with regards to chemical nature and relative elemental composition of the
13 samples. For example reference materials USGS40 and USGS41a are both glutamic acid but
14 of different ^{13}C and different ^{15}N abundance (Table 2). As amino acid the elemental C/N
15 ratio of glutamic acid resembles the C/N ratio typically associated with that of peptides and
16 proteins. Similarly, for 2H or ^{13}C analysis of hydrocarbons or long chain fatty acids reference
17 material pairs USGS67 / USGS69 or USGS70 / USGS72 respectively are ideally matched
18 and can thus serve as appropriate scale anchors (Table 2). For bulk 2H analysis of organic
19 compounds, in particular solids by TC/EA-IRMS (Thermal Conversion/Elemental Analyser -
20 IRMS) the need for well matched reference materials was a major concern up until recently.
21 Thanks to the efforts by Haiping Qi and Tyler Coplen (USGS Reston, VA, USA) scale
22 anchoring primary reference materials VSMOW and SLAP2 have become available sealed in
23 cold-welded silver tubes. This makes them thus amenable to be loaded on autosamplers for
24 solids and to be analysed side-by-side with solid samples for both bulk 2H and bulk ^{18}O
25 abundance by TC/EA-IRMS [21]. Recognising the resource constraints put on use of
26 VSMOW and SLAP2, secondary reference waters have been developed, also available sealed
27 in silver tubes from USGS Reston, whose δ^2H_{VSMOW} and $\delta^{18}O_{VSMOW}$ values cover the range
28 of the VSMOW/SLAP scale (Tables 2, 3).
29
30
31
32
33
34
35
36
37
38
39
40
41
42
43
44
45
46
47
48
49
50
51
52
53
54
55
56
57
58
59
60

5. Why appropriate choice of reference materials is crucial

The importance of which pair of reference materials are chosen to act as end-points, i.e. scale anchors cannot be stressed enough. The narrower the δ value range or bracket covered by one's choice of reference materials is, the less accurate the resulting normalisation will be of measured δ values outside that bracket (Figure 3). In the example shown in Figure 3 the $\delta^2\text{H}_{\text{VSMOW}}$ values of reference materials IAEA-CH-7 and IA-R002 bracket merely a $\delta^2\text{H}$ value range of 10.9 ‰, i.e. nowhere near the $\delta^2\text{H}$ value range of 427.5 ‰ covered by VSMOW / SLAP2. As a consequence, the linear correction equation based on IAEA-CH-7 and IR-R002 is significantly different from the scale normalisation equation obtained from contemporaneously analysed VSMOW and SLAP2 yielding as it does a "calibrated" $\delta^2\text{H}$ value of -198.0 ‰ for GISP that compared to its accepted $\delta^2\text{H}_{\text{VSMOW}}$ value of -189.7 ‰ is 8.3 ‰ too negative.

An analogous impact of too narrow a $\delta^{13}\text{C}$ value range is illustrated by the example shown in Figure 4a. Here the $\delta^{13}\text{C}$ range of 1.38 ‰ covered by reference materials IAEA-600 and USGS40 chosen as end points has a significant impact on accuracy of "corrected" $\delta^{13}\text{C}$ values. The difference of 0.283 ‰ may seem small between the "not-to-scale-calibrated" $\delta^{13}\text{C}$ value of -10.733 ‰ for IAEA-CH-6 and its scale normalised $\delta^{13}\text{C}$ value of -10.45 ‰. However, based on the acceptance criterion for modern IRMS instruments of 0.06 ‰ for repeatability of $\delta^{13}\text{C}$ measurements, the difference of 0.283 ‰ of (accepted minus corrected) is 4.72 times higher than this and thus statistically significant.

Considering how close proximity of two end-points results in pronounced differences for slope and off-set between the corresponding linear regression equation and that obtained from using appropriate scale anchors, one can understand why so-called single point calibrations cannot be and are not fit for purpose (Figure 4b). Obviously, a single point is not sufficient to define a linear equation given the slope s of a line is given by the ratio of

1
2
3 (y₂ - y₁)/(x₂ - x₁). A so-called single point "calibration" does most certainly not meet the
4
5
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22
23
24
25
26
27
28
29
30
31
32
33
34
35
36
37
38
39
40
41
42
43
44
45
46
47
48
49
50
51
52
53
54
55
56
57
58
59
60

(y₂ - y₁)/(x₂ - x₁). A so-called single point "calibration" does most certainly not meet the
aforementioned IUPAC requirement of normalising measured δ values to the size of the
relevant δ scale. Furthermore, by using merely one reference material as a single point of
comparison to correct measured δ values one either assumes a slope equal to 1 ($y = x + b$) or
an off-set of zero ($y = a \times x$). In case of the former the off-set "b" would be determined by
($RM_{\text{accepted}} - RM_{\text{measured}}$). In case of the latter the correction factor "a" would be determined
by ($RM_{\text{accepted}}/RM_{\text{measured}}$). However, either approach does not yield accurate, let alone scale
normalised δ values as Figure 4b illustrates for the scenario $y = x + b$. While the difference of
0.16 ‰ between single-point corrected and scale normalised $\delta^{13}\text{C}$ value of IAEA-CH-6 seems
small, it is still 2.67 times higher than the acceptance criterion of 0.06 ‰ for repeatability of
 $\delta^{13}\text{C}$ measurements.

The examples discussed above illustrate quite clearly the importance of ensuring the
two reference materials chosen as scale anchors are indeed end-points whose δ values extend
to if not match the size of the relevant δ scale. Furthermore, in addition to meeting with this
condition one should also take care when selecting reference materials with regards to their
physical properties and chemical elemental composition.

Modern on-line interfaces perform chemical conversion into a mixture of simple gases
and subsequent chromatographic separation of this gas mixture into individual gas peaks.
The combination of these different on-line processes in conjunction with transporting
permanent gases of low molecular weight in a carrier gas stream may be subject to potential
mass discriminatory effects. Further factors affecting measurement accuracy of δ values are
memory effects but also peak detection and peak integration by proprietary instrument
software. Adherence to the Identical Treatment principle of sample and standard can
minimise any analytical bias only when samples and reference materials are of comparable
chemical and elemental compositions, which is termed matrix matching.

1
2
3 For example, stable isotope analysis of ammonia salts or nitrogen bearing organic
4 compounds for their ^{15}N abundance is one example where appropriate choice of reference
5 materials in terms of matrix matching is extremely important. In ammonium salts or
6 organic amines or amides, nitrogen is present in states of oxidation ranging from -IV to -I.
7 However, reference materials USGS32 and USGS34 are both nitrates where nitrogen is
8 present in its highest state of oxidation, namely +V. While in ammonium salts and
9 nitrogen bearing organic compounds nitrogen will be oxidised to states of oxidation of 0
10 (N_2) or possibly +II (NO), nitrogen in nitrates cannot be oxidised any further but has to be
11 reduced to +IV (NO_2) or +II (NO). Therefore, organic reference materials such as L-
12 valine or caffeine should be chosen as scale anchors to ensure Identical Treatment
13 during compound conversion of nitrogen bearing organic compounds such as amino
14 acids or hetero-aromatic compounds respectively (Table 2).

15
16
17
18
19
20
21
22
23
24
25
26
27
28
29
30
31 Conversely, when analysing nitrate samples for their ^{15}N abundance, it is equally
32 important to use nitrate reference materials as scale anchors. This is not only important from
33 a matrix matching point of view but also from a sample conversion point of view. The
34 conversion of nitrates into ultimately N_2 gas relies on thermolysis of $[\text{NO}_3]^-$ into NO, NO_2
35 and O_2 , i.e. a reaction that generates oxygen rather than combustion, which consumes and
36 thus requires oxygen. This means conversion of nitrates in an EA has to be carried out
37 without the otherwise customary pulse of oxygen [22].

38
39
40
41
42
43
44
45
46
47 Last but not least attention is drawn to the perhaps confusing situation concerning
48 scale normalisation of $\delta^{18}\text{O}$ values. For historic reasons, there are three δ scales in use;
49 $\delta^{18}\text{O}_{\text{VSMOW}}$, $\delta^{18}\text{O}_{\text{VPDB}}$ and $\delta^{18}\text{O}_{\text{Air-O}_2}$. While reports of $\delta^{18}\text{O}_{\text{Air-O}_2}$ values are few and far
50 between, reporting $\delta^{18}\text{O}$ values of carbonates as $\delta^{18}\text{O}_{\text{VPDB}}$ values is still quite common even
51 though CIAAW's list of ^{18}O reference materials now states $\delta^{18}\text{O}_{\text{VSMOW}}$ values throughout,
52 including carbonate reference materials. It is crucially important to stay with one $\delta^{18}\text{O}$
53
54
55
56
57
58
59
60

reference scale when reporting and dealing with $\delta^{18}\text{O}$ values, and especially, when performing calculations involving $\delta^{18}\text{O}$ values. This may seem obvious, however, a glance at published articles reveals it is not [23, 24]. Subtracting $\delta^{18}\text{O}_{\text{VSMOW}}$ values from $\delta^{18}\text{O}_{\text{VPDB}}$ values makes about as much sense as subtracting kJ from kcal or ounces from grams. For meaningful scale normalisation of measured $\delta^{18}\text{O}$ values to either the VSMOW or the VPDB scale, accepted $\delta^{18}\text{O}$ values of the two reference materials chosen as scale anchors must be expressed on the same reference scale. Should instrumental or experimental set-up make it necessary for measurement results of two different sets of samples to be expressed on either of the two reference scales, this should be clearly stated in reports and publications. Furthermore, prior to any calculation involving both data sets $\delta^{18}\text{O}_{\text{VSMOW}}$ values have to be converted into $\delta^{18}\text{O}_{\text{VPDB}}$ values or *vice versa*. The conversion equations originally given by Friedman and O'Neil [25] were:

$$\delta^{18}\text{O}_{\text{VSMOW}} = 1.03086 \times \delta^{18}\text{O}_{\text{VPDB}} + 30.86 \quad (3)$$

$$\delta^{18}\text{O}_{\text{VPDB}} = 0.97006 \times \delta^{18}\text{O}_{\text{VSMOW}} - 29.94 \quad (4)$$

Please, note in equations (3) and (4) scale reference points SMOW and PDB in [25] have been amended to read VSMOW and VPDB in line with latest IUPAC guidelines [10].

However, with continuously advancing technology, new studies and resulting new insights, relations such as expressed in equations (3) and (4) are continuously being refined. Readers are therefore advised to visit e.g. the CIAAW website to check for the latest information [26]. With regards to the relation between $\delta^{18}\text{O}_{\text{VSMOW}}$ and $\delta^{18}\text{O}_{\text{VPDB}}$ values the CIAAW web page for oxygen reference materials currently states this equation [7,27]:

$$\delta^{18}\text{O}_{\text{VPDB}} = 0.97001 \times \delta^{18}\text{O}_{\text{VSMOW}} - 29.99 \quad (5)$$

6. Conclusions

- Most published methods sections describing stable isotopic analytical approaches and resulting δ values reflect shortcomings or even blatant misconceptions with regard to the proper use of reference materials, calibration, isotopic scale normalisation, and adherence to SI-rules and recommended IUPAC isotopic nomenclature.
- This paper offers detailed guidance about the mandatory use of isotopic reference materials for 2 end-point scale normalisation and calibrations against common isotopic scales.
- Practical examples of small data sets with various underlying scale normalisation and calibration strategies demonstrate the severity of isotopic artefacts resulting from inappropriately chosen "end"-points, inadmissible single-point calibration and violations of the principle of Identical Treatment of sample and standard.
- We offer a list of prominent suppliers of light stable isotope reference materials to best match sample matrices encountered by stable isotope practitioners with chemically similar reference materials to adhere to the Identical Treatment principle of sample and standard (Table 3).
- Proper use of reference materials and reporting of δ values fosters interdisciplinary communication while minimising ambiguity and reducing the jargon of scientific sub-disciplines.

7. References

- [1] de Groot PA, editor. Handbook of Stable Isotope Analytical Techniques. Amsterdam: Elsevier B.V.; 2004.
- [2] Fry B. Stable Isotope Ecology. New York (NY 10013): Springer; 2006.
<https://doi.org/10.1007/0-387-33745-8>; <http://www.jlakes.org/ch/web/0-387-33745-8.pdf>
- [3] Carter JF, Chesson LA, editors. Food Forensics: Stable Isotopes as a Guide to Authenticity and Origin. Boca Raton (FL): CRC Press, Taylor & Francis Group; 2017.
- [4] Sharp Z. Principles of Stable Isotope Geochemistry, 2nd Edition. University of New Mexico, Albuquerque (NM); 2017. http://digitalrepository.unm.edu/unm_oer/1/
- [5] Meier-Augenstein W. Stable Isotope Forensics - 2nd Edition: Methods and Forensic Applications of Stable Isotope Analysis. Chichester (UK): John Wiley & Sons Ltd; 2018.
- [6] Brand WA, Coplen TB. Stable isotope deltas: tiny, yet robust signatures in nature. *Isotopes Environ Health Stud.* 2012; 48: 393-409.
<https://doi.org/10.1080/10256016.2012.666977>
- [7] Brand WA, Coplen TB, Vogl J, et al. Assessment of international reference materials for isotope-ratio analysis (IUPAC Technical Report). *Pure Appl Chem.* 2014; 86: 425-467. <https://doi.org/10.1515/pac-2013-1023>
- [8] Coleman M, Meier-Augenstein W. Ignoring IUPAC guidelines for measurement and reporting of stable isotope abundance values affects us all. *Rapid Commun Mass Spectrom.* 2014; 28: 1953-1955. <https://doi.org/10.1002/rcm.6971>
- [9] Brand WA, Assonov SS, Coplen TB. Correction for the ^{17}O interference in $\delta(^{13}\text{C})$ measurements when analyzing CO_2 with stable isotope mass spectrometry (IUPAC Technical Report). *Pure Appl Chem.* 2010; 82: 1719-1733.
<https://doi.org/10.1351/PAC-REP-09-01-05>
- [10] Coplen TB. Guidelines and recommended terms for expression of stable-isotope-ratio and gas-ratio measurement results. *Rapid Commun Mass Spectrom.* 2011; 25: 2538-2560. <https://doi.org/10.1002/rcm.5129>
- [11] Park R, Epstein S. Metabolic fractionation of C^{13} & C^{12} in plants. *Plant Physiol.* 1961; 36: 133-138. <https://doi.org/10.1104/pp.36.2.133>

- 1
2
3 [12] Park R, Epstein S. Carbon isotope fractionation during photosynthesis. *Geochim*
4 *Cosmochim Acta*. 1960; 21: 110-126. [https://doi.org/10.1016-](https://doi.org/10.1016/S0016-7037(60)80006-3)
5 [7037\(60\)80006-3](https://doi.org/10.1016/S0016-7037(60)80006-3)
6
7
8 [13] Bender MM. Variations in the $^{13}\text{C}/^{12}\text{C}$ ratios of plants in relation to the pathway of
9 photosynthetic carbon dioxide fixation. *Phytochem*. 1971; 10: 1239-1244.
10 [https://doi.org/10.1016/S0031-9422\(00\)84324-1](https://doi.org/10.1016/S0031-9422(00)84324-1)
11
12 [14] Coplen TB. Normalization of oxygen and hydrogen isotope data. *Chem Geol*. 1988; 72:
13 293-297. [https://doi.org/10.1016/0168-9622\(88\)90042-5](https://doi.org/10.1016/0168-9622(88)90042-5)
14
15 [15] Horita J, Kendall C. Stable isotope analysis of water and aqueous solutions by
16 conventional dual-inlet mass spectrometry. In: de Groot PA, editor. *Handbook of Stable*
17 *Isotope Analytical Techniques, Volume I, Chapter 1*. Amsterdam
18 (Netherlands): Elsevier B.V.; 2004; pp. 1-37. [https://doi.org/10.1016/](https://doi.org/10.1016/B978-044451114-0/50003-X)
19 [B978-044451114-0/50003-X](https://doi.org/10.1016/B978-044451114-0/50003-X)
20
21 [16] Schimmelmann A, Qi H, Coplen TB, et al. New organic reference materials for
22 hydrogen, carbon, and nitrogen stable isotope-ratio measurements: caffeines, *n*-alkanes,
23 fatty acid methyl esters, glycines, L-valines, polyethylenes, and oils. *Anal Chem*. 2016;
24 88: 4294-4302. <https://doi.org/10.1021/acs.analchem.5b04392>
25
26 [17] Coplen TB. New guidelines for reporting stable hydrogen, carbon, and oxygen isotope-
27 ratio data. *Geochim Cosmochim Acta*. 1996; 60: 3359-3360.
28 [https://doi.org/10.1016/0016-7037\(96\)00263-3](https://doi.org/10.1016/0016-7037(96)00263-3)
29
30 [18] Coplen TB, Brand WA, Gehre M, et al. New guidelines for $\delta^{13}\text{C}$ measurements. *Anal*
31 *Chem*. 2006; 78: 2439-2441. <https://doi.org/10.1021/ac052027c>
32
33 [19] Werner RA, Brand WA. Referencing strategies and techniques in stable isotope ratio
34 analysis. *Rapid Commun Mass Spectrom*. 2001; 15: 501-519.
35 <https://doi.org/10.1002/rcm.258>
36
37 [20] Renpenning J, Kümmel S, Hitzfeld KL, et al. Compound-specific hydrogen isotope
38 analysis of heteroatom-bearing compounds via gas chromatography -chromium-based
39 High-Temperature Conversion (Cr/HTC) - Isotope Ratio Mass Spectrometry. *Anal Chem*.
40 2015; 87: 9443-9450.
41 <https://doi.org/10.1021/acs.analchem.5b02475>
42
43 [21] Qi H, Gröning M, Coplen TB, et al. Novel silver-tubing method for quantitative
44 introduction of water into high-temperature conversion systems for stable hydrogen and
45 oxygen isotopic measurements. *Rapid Commun Mass Spectrom*. 2010; 24: 1821-1827.
46 <https://doi.org/10.1002/rcm.4559>
47
48
49
50
51
52
53
54
55
56
57
58
59
60

- 1
2
3 [22] Gentile N, Rossi MJ, Delémont O, et al. $\delta^{15}\text{N}$ measurement of organic and
4 inorganic substances by EA-IRMS: a speciation-dependent procedure. *Anal Bioanal*
5 *Chem.* 2013; 405: 159-176. <https://doi.org/10.1007/s00216-012-6471-z>
6
7
8 [23] Dorval E, Piner K, Robertson L et al. Temperature record in the oxygen stable
9 isotopes of Pacific sardine otoliths: Experimental vs. wild stocks from the Southern
10 California Bight. *J Exp Mar Biol Ecol.* 2011; 397; 136-143.
11
12 <https://doi.org/10.1016/j.jembe.2010.11.024>
13
14 [24] Sakamoto T, Komatsu K, Yoneda M et al. Temperature dependence of $\delta^{18}\text{O}$ in
15 otolith of juvenile Japanese sardine: Laboratory rearing experiment with micro-
16 scale analysis. *Fish Res.* 2017; 194; 55-59.
17
18 <https://doi.org/10.1016/j.fishres.2017.05.004>
19
20 [25] Friedman I, O'Neil JR. Compilation of stable isotope fractionation factors of
21 geochemical interest. In: *Data of Geochemistry, Sixth Edition.* Geological Survey
22 Professional Paper 440-KK. Washington (DC): U.S. Government Printing Office;
23 1977; pp. KK1-KK12. <https://pubs.usgs.gov/pp/0440kk/report.pdf>
24
25 [26] <http://ciaaw.org/oxygen-references.htm>
26
27 [27] Kim S-T, Coplen TB, Horita J. Normalization of stable isotope data for carbonate
28 minerals: Implementation of IUPAC guidelines. *Geochim Cosmochim Acta.* 2015;
29 158; 276-289. <https://doi.org/10.1016/j.gca.2015.02.011>
30
31
32
33
34
35
36
37
38
39
40
41
42
43
44
45
46
47
48
49
50
51
52
53
54
55
56
57
58
59
60

1
2
3 Table 1. Difference between measured $\delta^2\text{H}$ values and $\delta^2\text{H}_{\text{VSMOW}}$ values calibrated by two
4 end-point scale normalization.
5
6
7

8
9 Table 2: Select list of reference materials that can serve as scale anchors (SCAN). Some
10 reference materials listed include a matching quality control (QC).^{a)}
11
12
13

14 Table 3. Distributors of light stable isotope reference materials (RMs) in alphabetical order.
15
16
17

18 Table 4. Generic batch sequence run sheet favouring high sample throughput under stable
19 experimental conditions using ^{13}C abundance analysis as example.
20
21
22
23
24
25
26
27

28
29 Figure 1. Comparison of variability and thus precision of quantitative mass spectrometric
30 analysis based (a) on absolute m/z abundance measurement and (b) on m/z abundance ratio
31 determination.
32
33
34
35
36
37

38 Figure 2. Flowchart for isotope abundance analysis observing the Identical Treatment
39 principle for both samples and reference materials serving as scale anchors.
40
41
42
43
44

45 Figure 3. Impact of reference material choice as end-points for $\delta^2\text{H}$ scale normalization on
46 correction equation and corrected $\delta^2\text{H}$ values.
47
48
49
50
51

52 Figure 4. Comparison of impact of reference material choice (a) as end-points or (b) as
53 single-point for $\delta^{13}\text{C}$ scale normalization on correction equation and corrected $\delta^{13}\text{C}$ values.
54
55
56
57
58
59
60

Tables for RM manuscript

Table 1: Difference between measured $\delta^2\text{H}$ values and $\delta^2\text{H}_{\text{VSMOW}}$ values calibrated by two end-point scale normalization.

Compound Name	$10^3 \times \delta^2\text{H}_{\text{measured}}$	$10^3 \times \delta^2\text{H}_{\text{VSMOW}}$	$10^3 \times \Delta \delta^2\text{H}$
VSMOW2	-0.36	0.00	0.36
NBS 22; $10^3 \times \delta^2\text{H}_{\text{accepted}} = -117.2$	-108.54	-117.00	8.46
SLAP2	-395.68	-427.50	31.82

Table 2: Select list of reference materials that can serve as scale anchors (SCAN). Some reference materials listed include a matching quality control (QC).^{a)}

Reference Material ID	Compound Name	$10^3 \times \delta^2\text{H}_{\text{VSMOW}}$	$10^3 \times \delta^{13}\text{C}_{\text{VPDB}}$	$10^3 \times \delta^{15}\text{N}_{\text{Air}}$
VSMOW2	Water (SCAN #1)	0.0		
GISP2	Water (QC1)	-258.3		
SLAP2	Water (SCAN #2)	-427.5		
USGS48	Water (SCAN #1)	-2.0 ± 0.4		
USGS47	Water (QC1)	-150.2 ± 0.5		
USGS49	Water (SCAN #2)	-394.7 ± 0.4		
USGS40	L-glutamic acid (SCAN #1)		-26.39 ± 0.04	-4.52 ± 0.06
USGS41a	L-glutamic acid (SCAN #2)		$+36.55 \pm 0.08$	$+47.55 \pm 0.15$
USGS64	Glycine (SCAN #1)		-40.81 ± 0.04	$+1.76 \pm 0.06$
USGS65	Glycine (QC)		-20.29 ± 0.04	$+20.68 \pm 0.06$
USGS66	Glycine (SCAN #2)		-0.67 ± 0.04	$+40.83 \pm 0.06$
USGS67	Hexadecane (SCAN #1)	-166.2 ± 1.0	-34.50 ± 0.05	
USGS68	Hexadecane (QC)	-10.2 ± 0.9	-10.55 ± 0.04	
USGS69	Hexadecane (SCAN #3)	$+381.4 \pm 3.5$	-0.57 ± 0.04	
USGS70	Icosanoic acid methyl ester (C20:0 FAME) (SCAN #1)	-183.9 ± 1.4	-30.53 ± 0.04	
USGS71	Icosanoic acid methyl ester (C20:0 FAME) (QC)	-4.9 ± 1.0	-10.50 ± 0.03	
USGS72	Icosanoic acid methyl ester (C20:0 FAME) (SCAN #2)	$+348.3 \pm 1.5$	-1.54 ± 0.03	
Reference Material ID	Compound Name		$10^3 \times \delta^{15}\text{N}_{\text{Air}}$	$10^3 \times \delta^{18}\text{O}_{\text{VSMOW}}$
IAEA-601	Benzoic acid (SCAN #1)			$+23.14 \pm 0.19$
IAEA-602	Benzoic acid (SCAN #2)			$+71.28 \pm 0.36$

1				
2				
3	USGS34	KNO_3 (^{15}N or ^{18}O SCAN #1)	-1.8 ± 0.1	-27.78 ± 0.37
4				
5	USGS32	KNO_3	$+180.0$	$+25.4 \pm 0.2$
6		(^{15}N SCAN #2; ^{18}O QC)		
7				
8	IAEA-NO-3	KNO_3 (^{18}O QC)		$+25.32 \pm 0.29$
9				
10	USGS35	NaNO_3	$+2.7 \pm 0.1$	$+56.81 \pm 0.31$
11		(^{15}N QC; ^{18}O SCAN #2)		
12				
13				

14 a) Reference material IDs and their δ values taken from: <http://www.ciaaw.org/reference-materials.htm>

15

16

17

18

19

20

21

22

23

24

25

26

27

28

29

30

31

32

33

34

35

36

37

38

39

40

41

42

43

44

45

46

47

48

49

50

51

52

53

54

55

56

57

58

59

60

Table 3: Distributors of light stable isotope reference materials (RMs) in alphabetical order.


Distributor	Web portal or contact	Comments	RM categories
AirLiquide (international)	http://isotope.airliquide-expertisecenter.com/	Gases like CO ₂ and SO ₂ from cylinders can only serve as monitoring gases in on-line applications	Air, trace gases in air, CO, CO ₂ , SO ₂ , NO _x , SF ₆ , hydrocarbons
ANSTO (Australia)	http://www.ansto.gov.au/ResearchHub/OurInfrastructure/A/CNS/CurrentResearch/ScientificHighlights/NDF-PE77/index.htm ; tde@ansto.gov.au	NDF-PE77 is isotopically indistinguishable from USGS77 powder	Polyethylene line NDF-PE77
Elemental Microanalysis (United Kingdom)	http://www.elementalmicroanalysis.com/product_list.php?top=IRMS%20supplement&category=204&sub=Certified	Website offers insufficient documentation of isotopic characterization (July 2018)	2 waters, 3 organic RMs
ERM [®] , European Reference Materials (Belgium)	https://crm.jrc.ec.europa.eu/c/By-application-field/Stable-isotopes/40476/	ERM [®] and IRMM RMs are identical and available from various vendors	Inorganic RMs
IAEA, International Atomic Energy Agency (Austria)	https://nucleus.iaea.org/rpst/referenceproducts/referencematerials/Stable_Isotopes/index.htm	Isotope data and inventory not up-to-date on website as of July 2018	Predominantly inorganic RMs, many waters
Indiana University, Department of Earth and Atmospheric Sciences (USA)	https://arndt.schimmelmann.us/welcome.html	e.g., Schimmelmann et al., 2016, <i>Analytical Chemistry</i> 88 , 4294-4302 http://dx.doi.org/10.1021/acs.analchem.5b04392	Organic RMs (gases, liquids, solids, GC-IRMS mixtures)
IRMM, Institute for Reference Materials and Measurements (Belgium)	https://crm.jrc.ec.europa.eu/c/By-application-field/Stable-isotopes/40476/	RMs from ERM [®] and IRMM are identical and available from various vendors	Inorganic RMs
Isometric Instruments (Canada)	http://www.isometricinstruments.com/gasstandards.html	Website offers insufficient documentation of isotopic characterization (July 2018)	Methane in air
NIST, National Institute of Standards and Technology (USA)	https://www-s.nist.gov/srmors/detail.cfm?ssearchstring=isotope	Isotope data and inventory not up-to-date on website as of July 2018	Predominantly inorganic RMs
NMI, National Measurement Institute (Australia)	chemref@measurement.gov.au ; available only for WADA-accredited forensic laboratories	Tobias and Brenna, 2018, <i>Drug Testing and Analysis</i> 10 (4), 781-785 https://doi.org/10.1002/dta.2309	Steroids for carbon stable isotope ratios only

OZTECH (USA)	isotopems@gmail.com http://www.si-science.co.jp/product/data/oztech.pdf	Website offers insufficient documentation of isotopic characterization (July 2018)	Pure gases CO ₂ , N ₂ , H ₂
Sercon Limited, Crewe (UK)	https://serconlimited.com/sercon_systems/standards/	Website offers insufficient documentation of isotopic characterization on RMs proprietary to Sercon (i.e. RMs with the prefix SC) (July 2018)	Inorganic and organic RMs; waters; flours and soils
SHOKO Science (Japan)	https://www.si-science.co.jp/global/en/index.html	RMs partially co-developed with JAMSTEC and Indiana University	Amino acids and waters
USGS, United States Geological Survey, Reston, Virginia (USA)	https://isotopes.usgs.gov/lab/referencematerials.html	Up-to-date isotope data on website; waters available in crimp-sealed silver capillary segments	Inorganic and organic RMs, including collagen and keratin
USGS, United States Geological Survey, Denver, Colorado (USA)	https://energy.usgs.gov/GeochemistryGeophysics/GeochemistryLaboratories/GasStandards.aspx	Dai et al., 2012, <i>Chemical Geology</i> 310-311 , 49-55. https://doi.org/10.1016/j.chemgeo.2012.03.008	Natural gas RMs of different geological origins


Table 4: Generic batch sequence run sheet favouring high sample throughput under stable experimental conditions using ^{13}C abundance analysis as example.

Sample ID ^{a)}	$10^3 \times \delta^{13}\text{C}$			$10^3 \times \delta^{13}\text{C}_{\text{VPDB}}$
	Measured	Corrected	Mean \pm S.D.	Accepted
BLANK				
RM1 (USGS40)				-26.39
RM1 (USGS40)				-26.39
RM1 (USGS40)				-26.39
RM1 (USGS40)				-26.39
<i>Samples k to l</i>				
AQC (IAEA-CH-6)				-10.45
AQC (IAEA-CH-6)				-10.45
AQC (IAEA-CH-6)				-10.45
AQC (IAEA-CH-6)				-10.45
<i>Samples m to n</i>				
RM2 (USGS41a)				+36.55
RM2 (USGS41a)				+36.55
RM2 (USGS41a)				+36.55
RM2 (USGS41a)				+36.55
BLANK				

a) RM = International Reference Material/s (here USGS40 and USGS41a) used to scale anchor measured $\delta^{13}\text{C}$ values on the VPDB scale. AQC = acquisition quality control sample/s (here IAEA-CH-6) used to control quality of scale normalisation.



1 Planning of analyses



Obtain ≥ 2 suitable international primary or secondary isotope reference materials (RMs) for your analytical application that should **bracket the δ -values that you expect** from your unknown samples. The **δ -value bracket size** should match the size of the δ -scale.

Optional: At a time when your stable isotope ratio analytical instruments work well, calibrate your own **tertiary laboratory RMs** against international RMs. The chosen RMs should bracket the δ -values that you expect from your unknown samples.

13 Performing analyses




Measure the ≥ 2 RMs and your unknown samples in a single analytical session in identical fashion, thus adhering to the **principle of Identical Treatment (IT) of standard and sample**. The resulting δ -values are raw, uncalibrated data.

20 Scale normalisation



Use proper statistical methods to **linearly scale-normalise raw δ -values** towards matching of δ -values of RMs with their prescribed δ -values.

27 Reporting of data




When describing analytical methods, explain which primary, secondary or tertiary RMs were used and **which δ -values were assigned to each RM**.


Propagate the errors of scale-normalised δ -values based on the empirically determined precision of repeat analyses and the accuracy of RM δ -values.

Use **SI-mandated and IUPAC-recommended nomenclature** when naming isotopes and expressing δ -values in terms of permil (‰) or milliurey (mUr).

1
2
3
4
5
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22
23
24
25
26
27
28
29
30
31
32
33
34
35
36



1
2
3
4
5
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22
23
24
25
26
27
28
29
30
31
32
33
34
35
36
37
38
39
40
41
42
43
44
45
46
47
48
49
50
51
52
53
54
55
56
57
58
59
60



a

IAEA-CH-6

-10

-20

-30


$\delta^{13}\text{C}$ corrected in ‰

$$\delta^{13}\text{C}_{\text{VPDB}} = 1.0109 \delta^{13}\text{C}_{\text{measured}} + 0.0963$$

$\delta^{13}\text{C}_{\text{VPDB}}$ of IAEA-CH-6 = -10.447 (accepted = -10.45 ‰)

$$\delta^{13}\text{C}_{\text{not-to-scale calib.}} = 0.9928 \delta^{13}\text{C}_{\text{measured}} - 0.3785$$

of IAEA-CH-6 = -10.733 (i.e. 0.283 ‰ too negative)



b

IAEA-CH-6

-10

-20

-30

-30

-20

-10

0

$\delta^{13}\text{C}$ measured in ‰

$$\delta^{13}\text{C}_{\text{VPDB}} = 1.0109 \delta^{13}\text{C}_{\text{measured}} + 0.0963$$

$\delta^{13}\text{C}_{\text{VPDB}}$ of IAEA-CH-6 = -10.447 (accepted = -10.45 ‰)

$$\delta^{13}\text{C}_{\text{corrected}} = \delta^{13}\text{C}_{\text{measured}} - 0.18$$

of IAEA-CH-6 = -10.61 (i.e. 0.16 ‰ too negative)